DCM064-5
Ed. 03/2012

fPSA

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BIBLIOGRAFIA

- Christensson A, Laurell CB, Lilja H., *Eur J Biochem*, 194, 755- 63 (1990).
- Watt KW, et. al., *Proc Nat Acad Sci USA*, 83, 3166-70 (1986).
- Chen Z, Prestigiacomo A, Stamey T., *Clin Chem*, 41, 1273-82 (1995).
- Wild D, *The Immunoassay Handbook*, Stockton Press (1994) p452.
- Junker R, Brandt B, Zechel C, Assmann G, *Clin Chem*, 43, 1588- 94 (1997).
- Prestigiacomo AF, Stamey TA, ' *Physiological variations of serum prostate antigen in the (4-10 ng/mL) range in male volunteers.*' *J. Urol* 1996;155:1977-80.
- Stamey TA, McNeal JE, Yemoto CM, Sigal BM, Johnstone IM, ' *Biological determinants of cancer progression in men with prostate cancer.*' *JAMA* 1999;281:1395-1400.
- Chen Z, Prestigiacomo A, Stamey T, ' *Purification and characterization of Prostate Specific Antigen (PSA) Complexed to a1-Anticymotrypsin: Potential reference Material for International Calibration of PSA Immunoassays.*' *Clin Chem*. 1995;41/9:1273-1282.
- Horton GL, Bahnsen RR, Datt M, Cfhan KM, Catalona WJ and Landenson JH. ' *Differences in values obtained with two assays of Prostate Specific Antigen.*' *J. Urol*. 1988;139:762-72.
- Stenman UH, Leinonen J, Alfthan H, Rannikko S, Tuhkanen Kand Alfthan O. ' *A complex between prostate specific antigen anda1-anticymotrypsin is the major form of prostate specific antigenin serum of patients with prostate cancer: assay of complex improves clinical sensitivity for cancer.*' *Cancer Res*. 1991;51:222-26.

Ed. 03/2012

DCM064-5

DiaMetra S.r.l. Headquater: Via Garibaldi, 18 –
20090 SEGRATE (MI) Italy
Tel. 0039-02-2139184 – 02-26921595
Fax 0039-02-2133354.

Manufactory: Via Pozzuolo 14, 06038 SPELLO (PG)
Italy
Tel. 0039-0742-24851
Fax 0039-0742-316197
E-mail: info@diametra.com

DCM064-5
Ed. 03/2012

fPSA

Direct immunoenzymatic determination of free Prostatic Specific Antigen (fPSA) in human serum or plasma

RUO

LOT

See external label

2°C  8°C

Σ = 96 tests

REF DKO064

INTENDED USE

fPSA kit is a direct solid phase enzyme immunoassay for quantitative determination of Free Prostatic Specific Antigen (fPSA) in human serum or plasma. fPSA kit is intended for laboratory use only.

1. CLINICAL SIGNIFICANCE

Prostate Specific antigen (PSA) is a serine protease with chymotrypsin-like activity. PSA is a single chain glycoprotein with a molecular weight of 28.4 kDA. PSA derives its name from the observation that it is a normal antigen of the prostate and is not found in any other normal or malignant tissue. PSA is released from the normal prostate and appears at low serum concentrations in healthy men. Studies with reverse transcription-PCR have shown that PSA also is expressed at a low concentration in peripheral blood cells and other tissues. High serum concentrations can be detected in patients with advanced prostate cancer (PCA). Therefore PSA is applied as a tumor marker for the clinical management of PCA. However, increased PSA concentrations in serum also occur in patients with benign prostate hyperplasia (BPH). Hence the goal in clinical laboratory is to discriminate clearly between BPH and PCA to spare the patient invasive diagnostic procedures, such as a prostate biopsy.

In human serum PSA occurs in two forms: free PSA (f-PSA) and complexed PSA. The major form is a complex of PSA and α_1 -antichymotrypsin (ACT). The fraction of f-PSA was shown to be substantially smaller in patients with untreated PCA than in patients with BPH. Therefore combined measurements of f-PSA and total PSA (t-PSA) may lead to a better discrimination between BPH and PCA. Some recent studies have already shown that the f-PSA/t-PSA ratio is helpful in the differential diagnosis of BPH and PCA.


PSA is found in benign, malignant and metastatic prostate cancer. Since prostate cancer is the second most prevalent form of male malignancy, the detection of elevated PSA levels plays an important role in the early diagnosis. Due to increased sensitivity, serum PSA levels have been found to be more useful than prostatic acid phosphatase (PAP) in the diagnosis and management of patients.

2. PRINCIPLE OF THE METHOD

Immunoenzymometric assay:

In this method, calibrators, patient specimens and/or controls (containing the native fPSA antigen) are first added to streptavidin coated wells. Biotinylated monoclonal and enzyme labeled antibodies are then added and the reactants mixed: these antibodies have high affinity and specificity and are directed against distinct and different epitopes of fPSA. Reaction between the various fPSA antibodies and native fPSA occurs in the microwells without competition or steric hindrance, forming a soluble sandwich complex

The interaction is illustrated by the following equation:

BtnAb(m) = Biotinylated Monoclonal Antibody (Excess Quantity)

AgfPSA = Native Antigen (Variable Quantity)


EnzAb = Enzyme labeled Antibody (Excess Quantity)

EnzAb -AgfPSA-BtnAb(m) = Antigen-Antibodies Sandwich Complex

k_a = Rate Constant of Association

k_{-a} = Rate Constant of Dissociation

Simultaneously, the complex is deposited to the well through the high affinity reaction of streptavidin and biotinylated antibody. This interaction is illustrated below:

StreptavidinC.W. = Streptavidin immobilized on well
Immobilized complex = sandwich complex bound to the well

After equilibrium is attained, the antibody-bound fraction is separated from unbound antigen by decantation or aspiration. The enzyme activity in the antibody-bound fraction is directly proportional to the native antigen concentration. The activity of the enzyme present on the surface of the well is

quantitated by reaction with a suitable substrate to produce colour. By utilizing several different calibrators of known antigen values, a dose response curve can be generated from which the antigen concentration of an unknown can be ascertained.

3. REAGENTS, MATERIALS AND INSTRUMENTATION

3.1. Reagents and materials supplied in the kit

1. fPSA Calibrators (6 vials, 1 mL each)

CAL0	REF DCE002/6406-0
CAL1	REF DCE002/6407-0
CAL2	REF DCE002/6408-0
CAL3	REF DCE002/6409-0
CAL4	REF DCE002/6410-0
CAL5	REF DCE002/6411-0
2. Conjugate (1 vial, 13 mL)
Antibodies anti fPSA conjugate with horseradish peroxidase (HRP) and anti fPSA biotinylated
REF DCE002/6402-0
3. Coated Microplate (1 microplate breakable)
Microplate coated with streptavidin
REF DCE002/6403-0
4. TMB Substrate (1 vial, 15 mL)
H₂O₂-TMB 0.26 g/L (*avoid any skin contact*)
REF DCE004-0
5. Stop Solution (1 vial, 15 mL)
Sulphuric acid 0.15 mol/L (*avoid any skin contact*)
REF DCE005-0
6. 50X Conc. Wash Solution (1 vial, 20 mL)
Phosphate buffer 50 mM pH 7.4; Tween-20 1 g/L
REF DCE006-0

3.2. Reagents necessary not supplied

Distilled water.

3.3. Auxiliary materials and instrumentation

Automatic dispenser.

Microplates reader (450 nm).

Notes

Store all reagents at 2-8°C in the dark.

Open the bag of reagent 3 (Coated Microplate) only when it is at room temperature and close it immediately after use; once opened, it is stable up to expiry date of the kit. Do not remove the adhesive sheets on the unutilised strips.

4. WARNINGS

- This kit is intended for in vitro use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- All human source material used in the preparation of the reagents has been tested and found negative for antibody to HIV 1&2, HbsAg, and HCV. No test method however can offer complete assurance that HIV, HBV, HCV or other infectious agents are absent. Therefore, the Calibrators

should be handled in the same manner as potentially infectious material.

- Some reagents contain small amounts of Proclin 300^R as preservative. Avoid the contact with skin or mucosa.
- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H₂O₂ to directed sunlight, metals or oxidants. Do not freeze the solution.
- This method allows the determination of fPSA from 0.5 ng/mL to 10.0 ng/mL.

5. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated. The reagents are stable until the expiry date when stored and handled as indicated.
- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately tested.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- Samples microbiologically contaminated, highly lipemic or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.

6. PROCEDURE

6.1. Preparation of the Calibrators (C₀...C₅)

The Calibrators are ready to use, were calibrated using a reference preparation (1st IS 96/670) and have the following concentration:

	C ₀	C ₁	C ₂	C ₃	C ₄	C ₅
ng/mL	0	0.5	1.0	2.5	5.0	10.0

Stability: until the expiry date printed on the label.

Once opened, the Calibrators are stable 6 months at 2-8°C.

6.2. Preparation of the Wash Solution

Dilute contents of wash buffer concentrate (50X) to 1000 mL with distilled or deionised water in a suitable storage container. For smaller volumes respect the dilution ratio of 1:50. The diluted buffer is stable at 2-8°C for at least 30 days.

6.3. Preparation of the Sample

- The specimens can be human serum or plasma; the usual precautions in the collection of venipuncture samples should be observed.
- For accurate comparison to established normal values, a fasting morning serum sample should be obtained.
- To obtain the serum, the blood should be collected in a venipuncture tube without additives or anti-coagulants. Allow the blood to clot. Centrifuge the specimen to separate serum from cells.
- Samples may be refrigerated at 2-8°C for a maximum period of 5 days. If specimens cannot be assayed within this time, the samples may be stored at temperatures of -20°C for up to 30 days.
- Avoid repetitive freezing and thawing.
- Patient specimens with fPSA concentrations above 10 ng/mL may be diluted (for example 1:10 or higher) with normal female serum (PSA = 0 ng/mL) and re-assayed. The sample's concentration is obtained by multiplying the result by the dilution factor.

6.4. Procedure

- **Allow all reagents to reach room temperature (22-28°C).**
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2-8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the calibration curve (C₀-C₅), two for each Control, two for each sample, one for Blank.

Reagent	Calibrator	Sample	Blank
Calibrator C ₀ -C ₅	50 µL		
Sample		50 µL	
Conjugate	100 µL	100 µL	
Incubate at room temperature (22-28°C) for 1 hour. Remove the contents from each well; wash the wells 3 times with 300 µL of diluted Wash Solution.			
TMB Substrate	100 µL	100 µL	100 µL
Incubate 15 minutes in the dark at room temperature (22±28°C).			
Stop Solution	100 µL	100 µL	100 µL
Shake the microplate gently. Read the absorbance (E) at 450 nm against Blank within 5 minutes.			

7. RESULTS

7.1. Validity of the assay

The OD of calibrator 5 should be ≥ 1.3.

7.2. Quality control

Each laboratory should assay controls at levels in the low, medium and high ranges of the dose response curve for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

7.3. OD conversion

The optical densities (O.D.s) of some calibrators and samples may be higher than 2.0; in such a case, they could be out of the measurement range of the microplate reader. It is therefore necessary, for O.D.s higher than 2.0, to perform a reading at 405 nm (= wavelength of peak shoulder) in addition to 450 nm (peak wavelength) and 620 (reference filter for the subtraction of interferences due to the plastic). For microplate readers unable to read the plate at 3 wavelengths at the same time, it is advisable to proceed as follows:

- Read the microplate at 450 nm and at 620 nm.
- Read again the plate at 405 nm and 620 nm.
- Find out the wells whose ODs at 450 nm are higher than 2.0
- Select the corresponding ODs read at 405 nm and multiply these values at 405 nm by the conversion factor 3.0 (where OD 450/OD 405 = 3.0), that is: $OD_{450\text{ nm}} = OD_{405\text{ nm}} \times 3.0$.

Warning: The conversion factor 3.0 is suggested only. For better accuracy, the user is advised to calculate the conversion factor specific for his own reader.

7.4. Data reduction - automated method

Use the 4 parameters logistic – preferred – or the smoothed cubic spline function as calculation algorithm.

NOTE: If computer controlled data reduction is used to calculate the results of the test, it is imperative that the predicted values for the calibrators fall within 10% of the assigned concentrations.

7.5. Data reduction: manual method

A dose response curve is used to ascertain the concentration of fPSA in unknown specimens.

- Record the OD obtained from the printout of the microplate reader as outlined in Example 1.
- Plot the OD for each duplicate calibrator versus the corresponding fPSA concentration in ng/mL on linear graph paper (do not average the duplicates of the calibrators before plotting).
- Draw the best-fit curve through the plotted points. To determine the concentration of fPSA for an unknown, locate the average absorbance of the duplicates for each unknown on the vertical axis of the graph, find the intersecting point on the curve, and read the concentration (in ng/mL) from the horizontal axis of the graph (the duplicates of the unknown may be averaged as indicated).

EXAMPLE 1				
Sample ID	Well number	OD	Mean OD	Value ng/mL
CAL 0	A1	0.019	0.019	0
	B1	0.021		
CAL 1	C1	0.118	0.120	0.5
	D1	0.121		
CAL 2	E1	0.211	0.214	1.0
	F1	0.217		
CAL 3	G1	0.703	0.701	2.5
	H1	0.699		
CAL 4	A2	1.617	2.338	5.0
	B2	1.661		
CAL 5	C2	2.393	2.338	10.0
	D2	2.282		
patient	E2	0.395	0.393	1.67
	F2	0.392		

The data presented in Example 1 are for illustration only and should not be used in lieu of a dose response curve prepared with each assay.

8. INTERPRETATION

- fPSA is elevated in benign prostatic hyperplasia (BPH). Clinically an elevated fPSA value alone is not of diagnostic value as a specific test for differential diagnosis of BPH. The ratio of fPSA/tPSA is a better marker and should be used in conjunction with other clinical observations (DRE) and diagnostic procedures (prostate biopsy).

- When the total PSA (tPSA) reads 4-10 ng/mL the fPSA/tPSA ratio is useful in the differential diagnosis of BPH and PC (Prostate Cancer). Depending on the ratio the probability can be determined as follows:

fPSA/tPSA Ratio	Probability of Prostate Cancer
0-10%	55%
10-15%	28%
15- 20%	25%
> 20%	10%

9. REFERENCE VALUES

Expected Values for the PSA Elisa Test System

PSA	
Healthy Males	≤ 1.3 ng/mL

Please pay attention to the fact that the determination of a range of expected values for a “normal” population in a given method is dependent on many factors, such as specificity and sensitivity of the method used and type of population under investigation. Therefore each laboratory should consider the range given by the Manufacturer as a general indication and produce their own range of expected values based on the indigenous population where the laboratory works.

Positive results should be verified concerning the entire clinical status of the patient. Also every decision for therapy should be taken individually. It is recommended that each laboratory establishes its own normal and pathological ranges of seric PSA.

10. PERFORMANCE AND CHARACTERISTICS

10.1. Precision

The within and between assay precision of the fPSA kit were determined by analyses on three different levels of control sera. The number, mean value, Calibrator deviation (SD) and coefficient of variation for each of these control sera are presented in Table 2 and Table 3.

Table 3 Between Assay Precision* (ng/mL)				
Sample	N	Mean	SD	C.V.
Level 1	10	0.52	0.04	11.3%
Level 2	10	2.34	0.22	5.8%
Level 3	10	7.70	0.68	5.2%

*As measured in ten experiments in duplicate

Table 2 Within Assay Precision (ng/mL)				
Sample	N	Mean	SD	C.V.
Level 1	20	0.43	0.04	9.3%
Level 2	20	2.57	0.20	7.8%
Level 3	20	8.20	0.73	8.9%

10.2. Accuracy-method comparison

Diametra fPSA kit was compared to a commercially available fPSA kit. 167 serum samples were tested.

The regression curve is:

Diametra = 0.9649*(reference kit) + 0.0189

R squared = 0.957

10.3. Sensitivity

fPSA kit has a sensitivity of 0.052 ng/mL.

10.4. Specificity

No interference was detected with the performance of Diametra PSA Elisa upon addition of massive amounts of the following substances to a human serum pool.

Compound	Concentration Added
AFP	10 µg/mL
Atropine	100 µg/mL
Acetylsalicylic Acid	100 µg/mL
Ascorbic Acid	100 µg/mL
Caffeine	100 µg/mL
Dexamethasone	10 µg/mL
Flutamide	100 µg/mL
hCG	100 IU/mL
hLH	100 IU/mL
Methotrexate	100 µg/mL
Prolactin	100 µg/mL
TSH	100 mIU/mL

11. WASTE MANAGEMENT

Reagents must be disposed off in accordance with local regulations.

BIBLIOGRAPHY

- Christensson A, Laurell CB, Lilja H., *Eur J Biochem*, 194, 755- 63 (1990).
- Watt KW, et. al., *Proc Nat Acad Sci USA*, 83, 3166-70 (1986).
- Chen Z, Prestigiacomo A, Stamey T., *Clin Chem*, 41, 1273-82 (1995).
- Wild D, *The Immunoassay Handbook*, Stockton Press (1994) p452.
- Junker R, Brandt B, Zechel C, Assmann G, *Clin Chem*, 43, 1588- 94 (1997).
- Prestigiacomo AF, Stamey TA, ' *Physiological variations of serum prostate antigen in the (4-10 ng/mL) range in male volunteers.*' *J. Urol* 1996;155:1977-80.
- Stamey TA, McNeal JE, Yemoto CM, Sigal BM, Johnstone IM,'*Biological determinants of cancer progression in men with prostate cancer.*' *JAMA* 1999;281:1395-1400.
- Chen Z, Prestigiacomo A, Stamey T, 'Purification and characterization of Prostate Specific Antigen (PSA) Complexed to a1-Anticymotrypsin:Potential reference Material for International Calibratorization of PSA Immunoassays.' *Clin Chem*. 1995;41/9:1273-1282.
- Horton GL, Bahnson RR, Datt M, Cfhan KM, Catalona WJ and Landenson JH.'*Differences in*

values obtained with two assays of Prostate Specific Antigen.' *J. Urol*. 1988;139:762-72.

- Stenman UH, Leinonen J, Alfthan H, Rannikko S, Tuhkanen Kand Alfthan O.'*A complex between prostate specific antigen anda1-anticymotrypsin is the major form of prostate specific antigenin serum of patients with prostate cancer:assay of complex improves clinical sensitivity for cancer.*' *Cancer Res*.1991;51:222-26.

Ed. 03/2012

DCM064-5

DiaMetra S.r.l. Headquater: Via Garibaldi, 18 – 20090 SEGRATE (MI) Italy

Tel. 0039-02-2139184 – 02-26921595

Fax 0039-02-2133354.

Manufactory: Via Pozzuolo 14, 06038 SPELLO (PG) Italy

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- Christensson A, Laurell CB, Lilja H., *Eur J Biochem*, 194, 755- 63 (1990).
- Watt KW, et. al., *Proc Nat Acad Sci USA*, 83, 3166-70 (1986).
- Chen Z, Prestigiacomo A, Stamey T., *Clin Chem*, 41, 1273-82 (1995).
- Wild D, *The Immunoassay Handbook*, Stockton Press (1994) p452.
- Junker R, Brandt B, Zechel C, Assmann G, *Clin Chem*, 43, 1588- 94 (1997).
- Prestigiacomo AF, Stamey TA, ' *Physiological variations of serum prostate antigen in the (4-10 ng/mL) range in male volunteers.*' *J. Urol* 1996;155:1977-80.
- Stamey TA, McNeal JE, Yemoto CM, Sigal BM, Johnstone IM, ' *Biological determinants of cancer progression in men with prostate cancer.*' *JAMA* 1999;281:1395-1400.
- Chen Z, Prestigiacomo A, Stamey T, ' *Purification and characterization of Prostate Specific Antigen (PSA) Complexed to a1- Anticymotrypsin: Potential reference Material for International Standardization of PSA Immunoassays.*' *Clin Chem*. 1995;41/9:1273-1282.
- Horton GL, Bahnson RR, Datt M, Cphan KM, Catalona WJ and Landenson JH. ' *Differences in values obtained with two assays of Prostate Specific Antigen.*' *J. Urol*. 1988;139:762-72.
- Stenman UH, Leinonen J, Alfthan H, Rannikko S, Tuhkanen Kand Alfthan O. ' *A complex between prostate specific antigen and a1-anticymotrypsin is the major form of prostate specific antigen in serum of patients with prostate cancer: assay of complex improves clinical sensitivity for cancer.*' *Cancer Res*. 1991;51:222-26.

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ERROR POSSIBLE CAUSES / SUGGESTIONS**No colorimetric reaction**

- no conjugate pipetted reaction after addition
- contamination of conjugates and/or of substrate
- errors in performing the assay procedure (e.g. accidental pipetting of reagents in a wrong sequence or from the wrong vial, etc.)

Too low reaction (too low ODs)

- incorrect conjugate (e.g. not from original kit)
- incubation time too short, incubation temperature too low

Too high reaction (too high ODs)

- incorrect conjugate (e.g. not from original kit)
- incubation time too long, incubation temperature too high
- water quality for wash buffer insufficient (low grade of deionization)
- insufficient washing (conjugates not properly removed)

Unexplainable outliers

- contamination of pipettes, tips or containers
- insufficient washing (conjugates not properly removed) too high within-run
- reagents and/or strips not pre-warmed to CV% Room Temperature prior to use
- plate washer is not washing correctly (suggestion: clean washer head)
- too high between-run - incubation conditions not constant (time, CV % temperature)
- controls and samples not dispensed at the same time (with the same intervals) (check pipetting order)
- person-related variation

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